# SoftCloud® Ace, SoftCloud® Pro

en	Alternating mattress systems User Manual .....	3
da	Alternierende madrassystem Brugsanvisning .....	21
es	Sistemas de colchón de aire alternante Manual del usuario .....	39
fi	Vaihtelevat patjajärjestelmät Käyttöohje .....	61
it	Materassi antidecubito a pressione alternata Manuale d'uso .....	79
no	Alternierende madrassystem Bruksanvisning .....	101
pt	Sistemas de colchão de pressão alternada Manual de utilização .....	119
sv	Madrassystem med växeltryck Bruksanvisning .....	141

This manual MUST be given to the user of the product.  
BEFORE using this product, this manual MUST be read and saved  
for future reference.

**Yes, you can.**

© 2023 Invacare Corporation

All rights reserved. Reproduction, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

Making Life's Experiences Possible is a registered trademark in the U.S.A.

---

# Contents

---

This manual MUST be given to the user of the product.  
BEFORE using this product, this manual MUST be read and saved  
for future reference.

<b>1 General</b>	<b>4</b>
1.1 Introduction	4
1.2 Symbols in this Document	4
1.3 Compliance	4
1.4 Warranty Information	4
1.5 Limitation of Liability	4
1.6 Service Life	4
<b>2 Safety</b>	<b>5</b>
2.1 Safety Information	5
2.2 Labels on the Product	6
2.3 Symbols on the product	6
<b>3 Product Overview</b>	<b>7</b>
3.1 Intended Use	7
3.1.1 Intended Users	7
3.1.2 Indications	7
3.1.3 Contraindications	7
3.2 Product Description	7
3.3 Control Unit	7
3.4 Mattress	7
<b>4 Setup</b>	<b>8</b>
4.1 Safety information	8
4.2 Installing the Mattress System	8
4.3 Activating the control unit	9
4.4 Removing the mattress system	9
<b>5 Usage</b>	<b>10</b>
5.1 Safety Information	10
5.2 Operating the Control Unit	10
5.3 Using the mattress system	11
5.4 Heel Cell Zones - SoftCloud Pro	11
5.5 CPR Procedure	11
5.6 Audible Warning Functions	12
5.7 Transporting a patient on the mattress	12
<b>6 Maintenance</b>	<b>13</b>
6.1 Inspection	13
6.2 Cleaning and Disinfection	13
6.2.1 General Safety Information	13
6.2.2 Cleaning Instructions	13
6.2.3 Disinfection Instructions	14
6.3 Service	14
6.4 Replacing the air filter	14
6.5 Replacing fuse	14
<b>7 After Use</b>	<b>15</b>
7.1 Storage	15
7.2 Reconditioning	15
7.3 Disposal	15
<b>8 Troubleshooting</b>	<b>16</b>
8.1 Identifying and Repairing Faults	16
<b>9 Technical Data</b>	<b>17</b>
9.1 Control Unit Specifications	17
9.2 Mattress specifications	17
9.3 Environmental Conditions	18
9.4 Guidance and Manufacturer's Declaration	18

# 1 General

## 1.1 Introduction

This user manual contains important information about the handling of the product. To ensure safety when using the product, read the user manual carefully and follow the safety instructions.

Only use this product if you have read and understood this manual. Seek additional advice from a healthcare professional who is familiar with your medical condition and clarify any questions regarding the correct use and necessary adjustment with the healthcare professional.

Note that there may be sections in this document, which are not relevant to your product, since this document applies to all available models (on the date of printing). If not otherwise stated, each section in this document refers to all models of the product.

The models and configurations available in your country can be found in the country-specific sales documents.

Invacare reserves the right to alter product specifications without further notice.

Before reading this document, make sure you have the latest version. You find the latest version as a PDF on the Invacare website.

If you find that the font size in the printed document is difficult to read, you can download the PDF version from the website. The PDF can then be scaled on screen to a font size that is more comfortable for you.

For more information about the product, for example product safety notices and product recalls, contact your Invacare distributor. See addresses at the end of this document.

In case of a serious incident with the product, you should inform the manufacturer and the competent authority in your country.

## 1.2 Symbols in this Document

Symbols and signal words are used in this document and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

### WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

### CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

### NOTICE

Indicates a hazardous situation that could result in damage to property if it is not avoided.

### Tips and Recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

## Other Symbols

(Not applicable for all manuals)

### UK Responsible Person

Indicates if a product is not manufactured in the UK.

### Triman

Indicates recycling and sorting rules (only relevant for France).

## 1.3 Compliance

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485.

This product features the CE mark, in compliance with the Medical Device Regulation 2017/745 Class I.

This product features the UKCA mark, in compliance with Part II UK MDR 2002 (as amended) Class I.

We are continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

We only use REACH compliant materials and components.

We comply with the current environmental legislations WEEE and RoHS.

## 1.4 Warranty Information

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business in the respective countries.

Warranty claims can only be made through the provider from whom the product was obtained.

## 1.5 Limitation of Liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorised modifications and/or use of unsuitable spare parts

## 1.6 Service Life

The expected service life of this product is five years when used daily and in accordance with the safety instructions, maintenance intervals and correct use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

## 2 Safety

### 2.1 Safety Information

Education, clinical judgement and action-based planning based on vulnerability are fundamental factors in prevention of pressure injuries.

A range of assessment scales can be used as a formal method of assessing risk from pressure injury development, and should be used in conjunction with an informal assessment (informed nursing judgement). Informal assessment is considered to be of greater importance and clinical value.

**WARNING!**

**Risk of Serious Injury or Damage**

Improper use of this product may cause injury or damage.

- If you are unable to understand the warnings, cautions or instructions, contact a health care professional or provider before attempting to use this equipment.
- Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manual, service manual or instruction sheet supplied with this product or optional equipment.

**WARNING!**

**Risk of fire or explosion!**

A cigarette can burn a hole in the bed surface and cause damage to the device. Also, patient clothing, bed sheets, etc. may be combustible and cause a fire. Failure to observe this warning can result in a severe fire, property damage and cause physical injury or death.

- Take special care in oxygen rich environments.
- Do not smoke.

**WARNING!**

**Risk of Developing Pressure Injuries**

The device delivers effective pressure redistribution, when the support surface is covered with a cotton, cotton combination or linen bed sheet, and any one of these is the only item deployed between the support surface and the user.

- Bed sheets must be loosely fitted, with creases smoothed out.
- Ensure that the support surface in contact with the user is kept free from crumbs and other food debris, and that drip cables, stents, and other foreign objects do not become entrapped between the support surface and the user.
- Heated over blankets must only be used in consultation with a qualified healthcare professional, as an increase in temperature can increase the risk of developing pressure injuries.

**CAUTION!**

**Risk of Injury**


Non-original or wrong parts may affect the function and safety of the product.

- Only use original parts for the product in use.
- Due to regional differences, refer to your local Invacare catalogue or website for available options or contact your Invacare distributor. See addresses at the end of this document.

## 2.2 Labels on the Product

The labels shown are for exemplification purposes alone. Note that the labels on your product might differ slightly from these examples.

### Control Unit

A	Product Name and Model
B	Input Voltage and Frequency Range
C	Max. Input Current
D	Fuse type
E	Ingress Protection (IP) Rating
F	Serial Number
G	Date of Manufacture
H	Unique Device Identification
I	Manufacturer
J	Country of Origin
K	European representative

### Mattress / Cover

Since the identification labels on the mattress and cover may change depending on model and date of manufacture, the labels for these components are not shown. For explanations of the symbols, which are printed on the mattress and cover, refer to section *Symbols on the product* in this document.

## 2.3 Symbols on the product

	European Conformity		UK Conformity Assessed
	Medical Device		European representative

	Manufacturer		Date of manufacture
	Unique Device Identification		Serial number
	Caution		CPR label
	Power		WEEE conform
	Type BF applied part		Class II equipment
	Functional Earth		User weight limit*
	Read the User Manual before using this product and follow all instructions for safety and use.		
	The background colour of this symbol is blue on product labels.		

	Do not pierce or cut		Do not iron
	Do not bleach		Tumble dry low temperature
	Do not put near flame		Line dry
	Hand wash		Do not dry clean
	Washing max. 95 °C (cover only)		

\* Minimum / maximum user weight per model, see chapter *Technical Data*.

## 3 Product Overview

### 3.1 Intended Use

This pressure redistribution mattress and control unit are intended to be used in conjunction with an appropriately sized bed frame.

The mattress when in dynamic mode offers alternating pressure support.

This product delivers effective pressure redistribution to users, when the support surface is covered with a cotton, cotton combination or linen bed sheet, and any one of these is the only item deployed between the support surface and the user.

#### 3.1.1 Intended Users

Intended users are adults and adolescents with restricted or limited mobility and or activity in a hospital, long term care or homecare setting.

#### 3.1.2 Indications

Suitable for supporting the management of all categories of pressure injuries when combined with an individual and comprehensive pressure injury protocol.

Suitable for use in home care, residential, nursing and acute care settings.


#### 3.1.3 Contraindications

Not suitable for users with unstable spinal cord fractures and/or cervical traction.

Always consult a physician prior to using this device.

### 3.2 Product Description

The following components are included within the scope of delivery:


(A)	Mattress System including mattress base, CPR tag, connector, handle with transport cap attached, top cover with zip.
(B)	Digital Control Unit
(C)	Carry Bag
(D)	User Manual
	Medical Grade Power Cord (not shown)

The control unit provides the air supply to the mattress. It is controlled via a touch panel. There is a visual and audible warning when pressure fails or power is interrupted.

The control unit redistributes the pressure by alternately inflating and deflating cells. This enables a frequent repositioning of the pressure.

### 3.3 Control Unit

- (A) Power button
- (B) Lock / Unlock button
- (C) Alternation failure indicator
- (D) Low pressure indicator
- (E) Service indicator
- (F) Power loss indicator
- (G) Audible warning mute button
- (H) Max. firmness button
- (I) Pressure settings
- (J) Mode button and indicator
- (K) Auto sensing function

### 3.4 Mattress

#### SoftCloud Ace

The mattress comprises air cells which all feature a permanently inflated internal cell to prevent the patient “bottoming out” in the event of low pressure.

The mattress includes three static head cells to provide static “pillow” support for optimum user comfort, while air pressure in the other cells is alternated.

The mattress also includes an independent heel zone of five micro cells for individualized support for this sensitive area.

#### SoftCloud Pro

The mattress comprises air cells which all feature a permanently inflated internal cell to prevent the patient “bottoming out” in the event of low pressure.

The mattress includes three static head cells to provide static “pillow” support for optimum user comfort, while air pressure in the other cells is alternated.

The mattress has a hinged back to allow the mattress to conform to the profile of the bed when articulated.

The mattress also includes an independent heel zone of five air cells for individualized support including permanent unloading for this sensitive area.

The mattress has side bolsters and the top cover has a double skirt with welded seams for maximum antibacterial protection.

## 4 Setup

### 4.1 Safety information

#### **WARNING!**

##### **Electrical shock hazard!**

- Do not remove control unit shroud.
- Refer to qualified service personnel.
- Before performing any maintenance to the control unit, disconnect the power lead from the wall outlet.
- Do not insert items into any openings of the control unit. Doing so may cause fire or electric shock by shorting the internal components.
- The control unit must be kept away from all heat sources and radiators during operation.
- Connect the equipment to a two- or three-prong wall outlet using the five meter power lead provided with the product.
- Position the device in such a way that ensures access to the power switch and CPR Connector at all times.

#### **WARNING!**

If this equipment is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.

- Do not modify this equipment without authorization of the manufacturer.

#### **WARNING!**

##### **Risk of entrapment!**

Patient entrapment with the bed side rails may cause injury or death. A thorough patient assessment should be completed and monitored and the equipment should be used as specified and maintained to reduce the risk of entrapment. Variations in bed side rail dimensions, and mattress thickness, size and density could increase the risk of entrapment.

- Mattress must fit bed frame and side rails to prevent patient entrapment. Follow the bed manufacturer's instructions.
- After any adjustments, repair or service and before use, make sure all attaching hardware is tightened securely. Side rails with dimensions different from the original equipment supplied or specified by the bed manufacturer may not be interchangeable and may result in entrapment or other injury.

#### **WARNING!**

This mattress is recommended to be installed on medical bed frames of an appropriate size with bed sides or assist rails, to accommodate the mattress. It is preferred that the rails are in the raised position whenever the patient is on the bed. Healthcare professionals assigned to each case should make the final decision whether assist rails are warranted after assessing patient risk of entrapment.

Controls on the foot bed end may be obstructed by the control unit on a few bed frames. It may be necessary to relocate the control unit.

- Before placing the patient on the bed, check that connector and power cord are clear of moving bed components.
- Operate all motorised bed frame functions through their full range of motion to be certain that there is no pulling, interference or pinching.
- Take care when positioning connector and cables to eliminate the risk of tripping hazards or strangulation.

#### **NOTICE!**

Inappropriate handling can damage the mattress cover.

- It is recommended that two people lift/carry mattresses.
- Avoid contact with jewellery, nails, abrasive surfaces etc.
- Do not drag mattresses.
- Avoid contact with wall, door frames, door catches or locks etc.
- Do not transport in roll cages unless completely protected from the sharp edges of the cage.

### 4.2 Installing the Mattress System

#### **NOTICE!**


- Remove all packaging before use.
- Keep all packing materials and instructions in the carry bag provided in case the product has to be shipped to an Invacare Service Centre.
- Carefully remove the control unit, mattress and accessories from the shipping cartons.
- Inspect all items for any damage that may have occurred during shipping. Any damaged or missing parts must be reported to an Invacare Service Centre immediately.


**NOTICE!**

- Ensure that the foot bed end is robust enough to support the control unit safely.
- Ensure there are no sharp objects which may come into contact with the mattress system.
- Check that the attachment of the mattress does not interfere with the movement or operation of the bed.
- Do not secure straps to bed side rails as straps will tear.
- Ensure that the positioning of the system does not interfere with the ability to disconnect the electrical power.
- Prior to activation ensure the CPR valve located near the head end of the mattress is fully engaged and that the mattress connector is firmly connected to the control unit.

1. Remove all covers, sheets and mattress from the bed.
2. Place the mattress directly on the frame of the bed.
3. Position the mattress top cover facing upwards and air hoses at foot of bed for control unit positioning.
4. On a standard bed, attach to the bed by securing the two adjustable straps under each end of the bed. Ensure buckles are securely fastened and straps are pulled tight. Or  
On a profiling bed, secure the two adjustable straps around the moveable sections of the bed.

**4.3 Activating the control unit**

1. Position the control unit by hanging the hooks over the foot board of the bed or side rails, ensuring that they are robust enough to hold the control unit.

**CAUTION!**

- Ensure the air hose does not kink between the bed frame and control unit.

2. Connect the mattress connector to the control unit.

**NOTICE!**

- Ensure that the power cord is in good condition.

3. Insert the power cord into the control unit then plug into a grounded 220V 50Hz electrical outlet.

4. Press the power button for approximately three seconds to activate the control unit.


- The pressure LEDs will flash indicating the system has activated.

**NOTICE!**

- Prior to pump activation ensure the CPR valve located near the head end of the mattress is fully engaged.
- Ensure that the mattress connector is firmly connected to the control unit.

5. Allow approximately 30 minutes for the mattress to fully inflate. Once ready, you should see the fourth pressure LED plus alternating mode LED illuminate to indicate that the system is ready for use (system automatically defaults to Alternating Mode after start-up). For further details see 5.3 *Using the mattress system, page 11*.
6. Once the Mattress is fully inflated the bedding can be placed. Fit sheets loosely enough to allow for free movement of the mattress air cells.
7. Ensure each sheet corner is placed through retaining buckle.

**4.4 Removing the mattress system**

1. Switch off the control unit and disconnect from mains supply (A).
2. Remove the mattress connector (B) from the control unit and disconnect the CPR tag (C).
3. Place control unit and power cord on top of the mattress and detach mattress from the bed frame.
4. Once air has been released from all cells, roll up the mattress and return all items to the carry bag for safe keeping.

**NOTICE!**

- Always follow the correct procedure when switching off the control unit. When the unit is switched off, it will issue a power loss alert. This alert can be canceled by pressing the power button for a minimum of two seconds.

**NOTICE!**

- Prior to restarting the system, ensure the CPR tag is replaced and all sealing connectors are firmly attached; and that the mattress connector is replaced and firmly connected to the control unit.
- Use target design to line up each plug with its corresponding socket.

## 5 Usage

### 5.1 Safety Information

#### WARNING!

Insufficient repositioning of the patient can lead to tissue compression and potential injury formation.

To relieve pressure, it is very important for the patients to reposition themselves, or to be repositioned, on a regular basis.

- Before using the product, always consult a qualified healthcare professional for clinical judgement.
- Monitor the patient frequently.

#### CAUTION!

- Make sure that the printed side of the mattress cover always faces upwards.
- Ensure that the distance between the surface of the mattress and the top of the side rail is at least 220 mm.

#### CAUTION!

##### Risk of damage to the mattress

If holes are present in the mattress cover, there is a risk that liquids may ingress and contamination may occur.

- Ensure that the mattress is not jammed or damaged by sharp edges.
- Do not place hypodermic needles, venflons, scalpels or other similarly sharp objects on or under the mattress.
- Do not use electrically heated blankets directly on or under the mattress.
- Ensure that all venflons are taped down correctly with no sharp edges exposed.
- When using bridging boards or other patient transfer aids, check for any sharp edges or burrs before use.
- When using the mattress on a profiling bed ensure that the knee break is used before the backrest.
- Attach medical equipment including infusion pumps and monitors to appropriate bed accessories.
- Avoid cigarette burns and pet claws puncturing the mattress cover.

### 5.2 Operating the Control Unit

#### Turning On and Off Power

1. Press power button for approximately three seconds to turn on or off the power.

#### Switching the Mode

1. Press mode button to switch between Alternation Mode (alternative cells cyclically inflating and deflating) and Static Mode (all cells fully inflated with no dynamic alternation). The bright LED indicates the active mode.
2. In order to activate or inactivate CLP function (Continuous Low Pressure), press mode button for two seconds.

#### Changing the pressure settings

Five pressure settings are available from soft (15 mmHg) to firm (55 mmHg) with steps of 10 mmHg.

When using the CLP function, the five pressure settings are from 15 mmHg to 23 mmHg with steps of 2 mmHg.

1. Press Plus button (+) to increase pressure.
2. Press Minus button (-) to decrease pressure.

The bright LEDs indicate the selected setting.

#### Locking and Unlocking buttons

1. Press Lock / Unlock button for approximately three seconds to lock all buttons on the control unit.  
A beep sounds and the amber LED illuminates to indicate the system is locked.
2. To unlock, press again for at least three seconds (beep sounds and amber LED turns off).

The control unit will automatically unlock in the event of a power failure.

#### Muting Audible Warning Signals

An audible warning signal can be muted while resolving the cause of the failure. If not resolved, the audible warning signal will resume after 17±3 minutes.

1. Press mute button to mute an audible warning signal.

#### Enabling Rapid Inflation to Maximum Pressure

1. Press Max. firmness button to enable rapid inflation to maximum pressure (60 mmHg).
2. To disable maximum pressure, press button again for three seconds.  
If not disabled, after 30±3 minutes the system will automatically revert back to alternation mode with medium pressure, for patient safety.

#### Auto Sensing (not applicable on all models)

The Auto sensing function is used when the mattress is transferred to a new user.

The system starts from a pressure reference value when calculating the individual pressure setting. When the user is lying on the mattress the differences between the “unloaded” and the “loaded” pressure are compared to a matrix and 1 of 5 pressure steps is selected. When adjusting the pressure the auto settings are overruled by the manual settings.

### 5.3 Using the mattress system

#### Establishing Pressure (supine / face up position)

1. The system is set up and inflated in static mode.
2. When the system is ready, the auto sensing function will indicate with a “beep” that the user can lay down on the mattress.
3. After 10-15 seconds the system has measured the impact the user have on the system and calculated a pressure setting.
4. The system is now ready to use and will turn to alternating therapy unless other selections are made.
5. If further adjustments are needed, the pressure settings can be altered with the pressure arrow buttons.
6. Assess whether the patient is comfortable and the system is functioning correctly by performing a “bottoming out” test.

 The system runs on the chosen setting but when the user changes and stays in a new position the pressure in the mattress can adjust (higher/lower). The sensor in the control unit will react and adjust the pressure if required to compensate for the new position.

#### “Bottoming Out” Test

When altering the pressure setting, ensure the patient is not “bottoming out” (insufficiently supported by the air cells and therefore coming in contact with bed base).

1. Ensure system is in alternation mode but is not undergoing an alternation.
2. With the patient lying in a supine position, unzip top cover just past sacral (bottom) region.
3. Slide your hand along a deflated cell under the patients sacral area (bottom). The inner static cell will remain inflated but your hand should slide easily between patient and base.
4. If hand can pass under patient then patient is adequately suspended and pressure can be lowered.
5. Repeat “bottoming out” test after pressure has been lowered.

 In the event of a system malfunction, the alarm will activate and pressure LEDs will flash.

#### Establishing Pressure (inclined position)

 The system runs on the chosen setting but when the user changes and stays in a new position the pressure in the mattress can adjust (higher/lower). The sensor in the control unit will react and adjust the pressure if required to compensate for the new position.

#### NOTICE!

- Return to the original pressure setting when the patient returns to the supine position, and perform a “bottoming out” test.
- Wait a minimum of 10-12 minutes between pressure adjustment and patient assessment, as it may take a full cycle for the system to adjust to any new setting.

### 5.4 Heel Cell Zones - SoftCloud Pro

For SoftCloud Pro, the five heel cells incorporate an option which allows the carer to permanently offload the heel of the user, ideal for the most vulnerable users such as those with limb ischaemia or existing wounds.

#### CAUTION!


##### Insufficient pressure relief


- Deflating more than one adjacent air cell can affect the support of the user during the normal alternating support cycle. This function should only be used for temporary procedures.

1. Open the push button flap **A** at the foot end of the mattress to access the 5 release connectors.
2. Disconnect the chosen connector **B** for permanent offloading.
  - 
    - Top connector is closest to the knee.
    - Connector 4–2 is below the knee and goes towards the heel.
    - Bottom connector is the last heel cell closest to the end of the mattress.
  -  It is not recommended to deflate more than one air cell directly beneath the area where relief is needed.
3. Reconnect the connector(s) to disengage permanent offloading. The air cell will be inflated again and return to normal operation mode.
4. Close the push button flap **A**.


### 5.5 CPR Procedure

1. Firmly pull the yellow rapid release CPR tag  from the side of the mattress to deflate the whole system.
2. Switch off the control unit.

 Mattress will start to deflate.

3.

When CPR is complete replace the yellow CPR Tag ensuring the sealing connectors are firmly attached and restart the control unit following chapter “Activating the control unit”. Using the target design to line up each plug with it's corresponding socket.

### ! NOTICE!

- Wait for the mattress system to gain optimal pressure.
- Perform a “bottoming out” test after inflating the mattress following rapid deflation.

## 5.6 Audible Warning Functions

The audible warning LED flashes, and an audible warning sounds, to indicate that the control unit or mattress pressure has failed. The LED will remain illuminated until appropriate pressure is restored. The audible warning can be silenced by pressing the audible warning mute button.

The system has four different audible warning signals and one service indicator, identified by five different symbols. The symbols are illustrated below.

Display	Audible Warning Signal	Description
	Alternating Mode failure	Indicates with an audible sound (press mute button to silence) and a flashing LED light when there is no alternation. The LED light will change to a fixed light within 5 min.
	Low Pressure	Indicates that pressure has fallen below the tolerance limit (below setting pressure –5 mmHg) or more. This will trigger the audible warning (press mute button to silence) and a flashing LED light. The LED light will change to a fixed light within 7 min.

Display	Audible Warning Signal	Description
	Power failure	Indicates no pressure output due to mains power failure. This will trigger the audible warning (press mute button to silence) and a flashing LED light. The LED light will change to a fixed light within 3 seconds. When the power is back on, the audible warning and the LED light will revert to normal mode.
	Initial failure (the mute button LED will be constantly lit)	Indicates that the system did not complete initialization within 50 minutes, the mute button LED light will have a fixed light.
	Service indicator	Indicates that service is needed, the LED light will have a fixed light. This will occur when the control unit has been running for approximately 9000 hours.

 If audible warning activates and the system fails to inflate or loses pressure, refer to section “Identifying and repairing faults”.

## 5.7 Transporting a patient on the mattress

1. Before patient transport, press the Mode button and wait at least 12 minutes for cells to inflate to maximum pressure.
  -  If a quicker response is required (up to 5-10 minutes), then use the Max Firmness button.
2. Once mattress pressure has reached maximum inflation, press the Power button to switch off the control unit. Switch off mains supply and unplug the power cord.
3. Remove the connector from the control unit and allow air to escape for a few seconds before inserting the transport cap into the connector air outlets to seal the system. This release softens the mattress surface for pressure relief and comfort. Ensure the cap is firmly secured across both air hose outlets.

### CAUTION!

**Air pressure is released from all internal static cells as well as alternating sections.**

- Regularly perform a “bottoming out” test to ensure the patient is appropriately supported.

## 6 Maintenance

### 6.1 Inspection

It is recommended to check mattresses (air cells and cover) for strike-through (this may include fluid ingress, stains, rips or damage) after the release of each patient or after each period of use by a suitably qualified and competent person.

#### Check mattresses

1. Unzip the cover completely.
2. Check for any staining on the white underside of the cover.
3. Check for any staining on the internal elements.
4. Replace any stained items and dispose of as per local authority procedure.

### 6.2 Cleaning and Disinfection

#### 6.2.1 General Safety Information

##### CAUTION!

##### Risk of Contamination

- Take precautions for yourself and use appropriate protective equipment.

##### CAUTION!

##### Risk of Electric Shock and Product Damage

- Switch off the device and disconnect from mains, if applicable.
- When cleaning electronic components consider their protection class regarding water ingress.
- Make sure that no water splashes to the plug or the wall outlet.
- Do not touch the power socket with wet hands.

##### NOTICE!

- Wrong fluids or methods can harm or damage the product.
- All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.
  - Never use corrosive fluids (alkalines, acid etc.) or abrasive cleaning agents. We recommend an ordinary household cleaning agent such as dishwashing liquid, if not specified otherwise in the cleaning instructions.
  - Never use a solvent (cellulose thinner, acetone etc.) that changes the structure of the plastic or dissolves the attached labels.
  - Always make sure that the product is completely dried before taking into use again.

For cleaning and disinfection in clinical or long-term care environments, follow your in-house procedures.

### Cleaning Intervals

##### NOTICE!

Regular cleaning and disinfection enhance smooth operation, increases the service life and prevents contamination.

Clean and disinfect the product

- regularly while it is in use,
- before and after any service procedure,
- when it has been in contact with any body fluids,
- before using it for a new user.

##### NOTICE!

- Keep a cleaning record as part of cleaning the system.

#### 6.2.2 Cleaning Instructions

##### NOTICE!

- The product does not tolerate cleaning in automatic washing plants, with high-pressure cleaning equipment or steam.

A cleaning record should be kept as part of cleaning the system.

#### Cleaning control unit

1. Wipe down the control unit casing and hose fittings with a damp cloth and suitable detergent.
2. Using a nylon brush, gently clean all crevices as they can harbour microorganisms.
3. Wipe off the control unit casing and hose fittings with a damp cloth to remove all detergent.
4. Air dry all treated surfaces.

#### Cleaning hangers

1. Wipe down the exterior of the hangers periodically, using a damp cloth and suitable detergent.

#### Cleaning covers

(Removal of contaminants such as dust and organic matter)

1. Remove all covers for laundering.
2. Launder the covers with maximum temperature as stated on product label, using a diluted detergent solution (Instructions on label).

##### NOTICE!

- Washing at higher temperatures will cause shrinkage.

#### Drying covers

1. Hang covers from a line or bar and drip dry in a clean indoor environment or tumble dry on a low heat setting.

##### NOTICE!

- Tumble dry setting must not exceed 40 °C.
- Do not tumble dry for longer than 10 minutes.
- Dry thoroughly before use.

### 6.2.3 Disinfection Instructions

- !** **NOTICE!**
- Only use disinfectants and methods approved by your local infection control institution and follow your local infection control policy.
  - Additional information on recommended disinfectants can be found on <https://vah-online.de/en/for-users>.

#### Disinfecting the Control Unit

1. Wipe down all generally accessible surfaces with a cloth and disinfectant.
2. Allow the product to air-dry.

#### Disinfecting the Air Cells

- !** **NOTICE!**
- Do not disconnect the pipes from individual air cells.
  - Do not machine wash or dry the air cells.

1. Disconnect air cells from the base by unfastening the press studs at each end.
2. Disconnect air pipes from main air hoses.
3. Slide each cell out from the cell straps.
4. Swab with a cloth damped with warm water containing an appropriate disinfectant (dependent on infection control protocols and local market requirements).
5. Dry thoroughly with a soft cloth before refastening.

### 6.3 Service

- !** **NOTICE!**
- The service indicator will be activated when the control unit has been running for approximately 9.000 hours, indicating that service is needed.

- !** **NOTICE!**
- Do not attempt to repair or service the control unit yourself! Always contact an authorised service technician.

### 6.4 Replacing the air filter

- !** **NOTICE!**
- Good filter maintenance is critical to maintain your mattress system in optimal operating condition. Failure to keep the filters clean will result in system downtime and increase repair costs. It is recommended that the air filter is replaced annually. Replacement air filters are available from an Invacare Service Centre.

1. Turn off the power supply to the control unit.
2. Disconnect the power cord and air hoses.
3. Place the control unit on a flat surface with its back panel uppermost (place soft cloth under unit to prevent scratches).
4. Carefully remove the air filter shroud. Remove and discard the filter and attach with new filter.
5. Re-attach the air filter shroud to the control unit.

The control unit is now ready for re-connection.

### 6.5 Replacing fuse

- CAUTION!**
- Ensure the replacement of fuses is carried out in accordance with local legislation.
  - Always use the same fuse rating when replacing a defect fuse.

1. Turn off the power supply to the control unit.
2. Remove the power cord from the electrical socket on the side of the control unit.
3. Insert a flat screwdriver into the small cutout above the power cord inlet.
4. Pull the screwdriver carefully towards to open the plastic carrier.
5. Remove and discard the “blown” fuse.
6. Insert a new fuse and close the plastic carrier.
7. Reverse steps four to one in order to reinstate the fuse carrier.

See also section “Activating the control unit” for more information.

## 7 After Use

### 7.1 Storage

#### NOTICE!

- Store device in a dry environment.
- Store device within a protective cover.
- Carefully roll up the device and store it on a clean, dry (off-floor) surface with no sharp edges to prevent possible damage.
- Never store other items on top of a device.
- Do not store device next to radiators or other heating devices.
- Protect device from direct sunlight.

Environmental conditions for storage, see “*Environmental conditions*”.

### 7.2 Reconditioning

This product is suitable for reuse. To recondition the product for a new user, carry out the following actions:

- Inspection according to service plan
- Cleaning and disinfection

For detailed information about inspection, cleaning and disinfection, refer to *6 Maintenance, page 13*.

### 7.3 Disposal

#### WARNING!

##### Environmental Hazard

Device contains batteries.

This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

- DO NOT dispose of batteries in normal household waste.
- DO NOT throw batteries into a fire.
- Batteries **MUST** be taken to a proper disposal site. The return is required by law and free of charge.
- Only dispose of discharged batteries.
- Cover terminals of lithium batteries prior to disposal.
- For information on the battery type see battery label or chapter *9 Technical Data, page 17*.

Be environmentally responsible and recycle this product through your recycling facility at its end of life.

Disassemble the product and its components, so the different materials can be separated and recycled individually.

The disposal and recycling of used products and packaging must comply with the laws and regulations for waste handling in each country. Contact your local waste management company for information.

## 8 Troubleshooting

### 8.1 Identifying and Repairing Faults

**WARNING!**
**Electric shock hazard!**

Opening the control unit could cause personal injury or equipment damage.  
– Do not try to open the control unit.

**NOTICE!**

– If the problem is not resolved, contact an Invacare provider or Invacare representative.

Warning/Fault	Cause	Solution
Control unit does not operate; no display lights	The Control Unit may not be attached to a power source  A fuse may need replacing in the control unit	<ol style="list-style-type: none"> <li>1. Check that the control unit is connected to mains power outlet with the correct voltage.</li> <li>2. Check that the control unit is switched on.</li> <li>3. Check the mains plug fuse. Check both control unit fuses (1 amp slow blow fuse) – fuses can be released using a screwdriver to push and turn.</li> </ol>
Patient is sinking or “bottoming out” whilst lying flat on the Mattress	The pressure may be set too low for the patient’s weight	<ol style="list-style-type: none"> <li>1. Increase the pressure setting by pressing up the pressure arrow.</li> <li>2. To check effective system performance, perform a “bottoming out” test. See section “Using the mattress system”.</li> </ol>
Audible Warning Mute LED lights up permanently  	Initial failure	<ol style="list-style-type: none"> <li>1. Mute the audible warning and turn off power.</li> <li>2. Check that the handle is intact, ensuring all four sealing connectors are firmly attached to control unit and air hoses. Check the CPR Tag is attached and all four sealing connectors are firmly secured.</li> <li>3. Check all air hoses along the inside of the mattress – each should be firmly connected. Check each air cell is securely attached to its connecting air pipe.</li> <li>4. Check all cells, pipes and hoses for any air leakage.</li> <li>5. Turn on Power.</li> </ol>
Low pressure indicator LED lights up    + audible warning	Pressure too low	<ol style="list-style-type: none"> <li>1. Mute the audible warning and turn off power.</li> <li>2. Check that the handle is intact, ensuring all four sealing connectors are firmly attached to control unit and air hoses. Check the CPR Tag is attached and all four sealing connectors are firmly secured.</li> <li>3. Check all air hoses along the inside of the mattress – each should be firmly connected. Check each air cell is securely attached to its connecting air pipe.</li> <li>4. Check all cells, pipes and hoses for any air leakage.</li> <li>5. Check that the air filter shroud is correctly secured and the air filter is clean.</li> <li>6. Turn on Power.</li> </ol>
Alternation Failure indicator LED lights up    + audible warning	Alternating Mode Failure (no alternation)	<ol style="list-style-type: none"> <li>1. Mute the audible warning and turn off power.</li> <li>2. Disconnect the air hoses to reduce pressure – reconnect when pressure has decreased.</li> <li>3. Restart the system by switching the power button off and on. The system automatically starts in Static Mode and will switch to Alternating Mode when operational pressure is reached.</li> <li>4. Reset the personal settings.</li> </ol>
Power loss indicator LED lights up    + audible warning	Power failure	<ol style="list-style-type: none"> <li>1. Mute the audible warning.</li> <li>2. Check that the power cable is firmly plugged into the mains outlet and the control unit. Check that the mains power is switched on.</li> <li>3. Check both control unit fuses (1 amp slow blow fuse) – fuse can be released using a screwdriver to push and turn.</li> </ol> <p>If power is restored within 20 minutes of failure, the system will run an initial start up sequence before returning to the last settings.</p>
Service Indicator LED lights up  	Service needed	<ol style="list-style-type: none"> <li>1. Contact an Invacare provider or Invacare representative.</li> </ol> <p> There is no risk for continued use even if the Service Indicator LED lights up.</p>

## 9 Technical Data

### 9.1 Control Unit Specifications

	Control Unit 1 in 3
Cycle Control	Purpose designed distributor valve supplying operating air to the inflatable cells
Cycle Time	12 minutes
Supply Voltage	220-240V 50/60Hz
Air Capacity	10 lpm
Noise level	Below 30dB
Fuse Rating	1A / 250V
Battery Source	VARTA, V80H, 1.2 V DC, 70 mAh
Power Rating	max. 1A
Control Unit Dimensions:	
Length	123 (±10) mm
Width	243 (±10) mm
Height	250 (±10) mm
Control unit weight	4.2 kg
Hose Connection	Push on connection handle
Emergency	CPR Tag
Mode of Operation	Non-continuous
Degree of protection	<b>IP21</b> — Protected against solid objects over 12.5 mm e.g. hands, large tools. Protected against vertically falling drops of water or condensation.

### 9.2 Mattress specifications

	SoftCloud Ace				SoftCloud Pro	
	Narrow	Wide	Large	XL	Narrow	Wide
Number of Cells	21 (cell in cell) including 3 static head cells and 5 heel cells				22 (cell in cell) including 3 static head cells and 5 heel cells	
Cell Height	200 mm				200 mm	
Dimensions:						
Length	2000 mm	2000 mm	2000 mm	2000 mm	2000 mm	2000 mm
Width	830 mm	880 mm	1030 mm	1160 mm	830 mm	880 mm
Height	200 mm	200 mm	200 mm	200 mm	200 mm	200 mm
Alternation type	1 in 3 cell cycles				1 in 3 cell cycles	
Mattress weight	8.0 kg	10.0 kg	11.2 kg	11.2 kg	9.8 kg	10.0 kg
Minimum/ maximum User Weight	0 – 220 kg		0 – 350 kg		0 – 220 kg	
Cell Material	0.15 mm TPU film laminated on 210 denier nylon fabric					
Base Material	Nylon fabric 420 denier with a 0.1 mm TPU coating					
Cover Material	100% Polyurethane surface, 100% Polyester inside					

### 9.3 Environmental Conditions

	Operation	Storage and Transport
Ambient temperature	10 °C – 40 °C	-10 °C – 60 °C
Relative humidity	30% – 75%, non-condensing	10% – 70%, non-condensing
Atmospheric pressure	70 – 106 kPa	70 – 106 kPa

### 9.4 Guidance and Manufacturer's Declaration

The mattress system is intended for use in the electromagnetic environment specified below. The customer or the user of the mattress system should ensure that it is used in such an environment.

#### Electromagnetic Emissions

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The mattress system use RF energy only for their internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The mattress system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

#### Electromagnetic Immunity

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s)	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	<5% $U_T$ (>95% dip in $U_T$ ) for 0.5 cycle 40% $U_T$ (60% dip in $U_T$ ) for 5 cycles 70% $U_T$ (30% dip in $U_T$ ) for 25 cycles <5% $U_T$ (>95% dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the mattress system requires continued operation during power mains interruptions, it is recommended that the mattress system be powered from an uninterruptible power supply or a battery.
Power frequency (50 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE:  $U_T$  is the a.c. mains voltage prior to application of the test level.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V <sub>rms</sub> 150 kHz to 80 MHz	3 V <sub>rms</sub>	Portable and mobile RF communications equipment should be used no closer to any part of the CT515, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  Recommended separation distance: $d = 1.167\sqrt{P}$ $d = 1.167\sqrt{P}$ 80 MHz to 800 MHz $d = 2.333\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).  Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a)</sup> , should be less than the compliance level in each frequency range <sup>b)</sup> .  Interference may occur in the vicinity of equipment marked with the following symbol:  
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>a)</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the mattress system is used exceeds the applicable RF compliance level above, the mattress system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the mattress system.

<sup>b)</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distances between portable and mobile RF communications equipment and the mattress system Alternating Control Unit**

The mattress system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the mattress system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the mattress system as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.167\sqrt{P}$	80 MHz to 800 MHz $d = 1.167\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.333\sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.689	3.689	7.379
100	11.667	11.667	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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Invacare New Zealand Ltd
Invacare AB
Invacare UK Operations Limited