

Combination Agreement



This document provides an overview of the verified compatibility between Invacare manual wheelchairs and selected third-party devices. It aims to help our customers in identifying combinations that are compliant with the EU Medical Device Regulation (MDR) 2017/745. All listed combinations have been assessed to help ensure safety, performance, and regulatory compliance when the products are used together as intended. This document is a reference guide only and may be updated periodically without notice as new assessments are completed. Please check with your Invacare sales representative or visit us at <https://www.invacare.eu.com/help-support> to make sure you have the latest version. The use of an Invacare manual wheelchair in combination with a third-party device not included in this document should be avoided until the responsible provider has conducted an independent risk assessment to verify compatibility between the wheelchair and such device.

● Agreement compliant with MDR2017/745

Invacare manual wheelchair	Including versions	Alber											Batec Mobility				Rehasense				Roam						Decon	Klaxon Mobility			Loopwheels LTD	Sunrise					Benoit Systemes													
		E-fix E35/E36	E-Motion M25	Twion T24	Viamobil V25	Viamobil eco V14	ViaGo	ViaPlus V12	Scalamobil S45	Scalamobil S38	Smoov Q10	e-pilot P15	Batec Mini 2	Batec Electric 2	Batec Rapid 2	Batec Hibrid 2	Paws City 12"/14"	Paws Cruiser 16"	Track Wheel C2 (double arm) carbon frame	Track Wheel A2 (double arm) aluminium frame	Eco Assist	Runner	E-Bike	Racer	Connect Carrier	Connect Aankop-pelstang	Connect Offroad	Adventus	Klaxon Klick Power (all versions)	Klaxon Klick Race (all versions)	Klaxon Klick Hybrid (all versions)	Loopwheels 2 (Urban, Extreme, LT)	Empulse Attitude (Junior / Manual / Hybrid / Power)	Empulse F35	Empulse F55	Empulse StreetJet 12"	Empulse WheelDrive	Empulse R20	Light Drive 3 / Minotor 3	Light Drive 3 XL / Minotor XL 3	Light Assist 3 / Pousseur 3	Light Assist 3 XL / Pousseur XL 3								
ACTION 1R		●	●	●	●	●	-	-	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-							
ACTION 2NG		●	●	●	●	●	●	-	●	●	-	-	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	●	●	●	-	-	-	-					
ACTION 2 TRANSIT LITE		-	-	-	●	●	●	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-							
ACTION 3JUNIOR	Evolutive	●	●	-	-	●	-	-	●	-	-	-	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	-	-					
	Light	●	●	●	●	●	-	-	●	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-					
ACTION 3NG	Comfort	●	●	●	●	●	●	-	●	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-					
	XT	●	●	●	●	●	●	-	●	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-					
ACTION 3NG: ROCKING CHAIR		-	-	-	-	-	-	-	-	-	●	-	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-				
ACTION 4NG	Comfort	●	●	●	●	●	●	-	-	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-				
ACTION 4NG: HEAVY DUTY VERSION	XLT	●	●		●	●	-	-	●	●	-	-	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-				
ACTION 5/ MYON HC		●	●	●	●	●	-	-	●	-	-	●	●	●	●	●	●*	●*	●*	●*	●*	●*	●	●	●	●	●	●	-	-	-	-	-	-	-	●*	●*	●	●	●	●	●	●	●	●	●				
ACTION 5/ MYON HC: RIGID VERSION		●	●	●	-	-	-	-	-	-	●	●	-	-	-	-	●*	●*	●*	●*	●*	●*	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-			
ACTION AMPLA		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-				
ESPRIT ACTION 4NG		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-				
ESPRIT ACTION		-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-				
SPIREA 4 NG		●	●	●	●	●	●	-	●	●	●	●	-	-	-	-	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	●	●	●			
REA AZALEA	Standard	●	●	●	●	-	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	●	●	●	●				
	Assist	●	●	●	●	-	-	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	-	-	-	-	-				
	Minor	●	●	●	●	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	●	-	-	-	-				
	Tall	●	●	●	●	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	●	-	-	-	-				
REA AZALEA BASE		●	●	●	●	-	-	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	●	●	●	●				
REA AZALEA MAX		●	●	●	●	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	●	●	●	●				
REA CLEMATIS		●	●	●	●	-	-	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	●	●			
REA CLEMATIS PRO	E-tilt	●	●	●	●	-	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	●	●			
REA DAHLIA	30, 45	●	●	●	●	-	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	●	●			
KÜSCHALL CHAMPION 2.0		-	●	●	-	●	-	-	-	-	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	●	●	●	-	●	●	●	●	-	●	●	●	●	●	●	●	●		
KÜSCHALL COMPACT 2.0		●	●	●	●	●	●	-	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-	●	●	●	-	-	-	-	●	●	-	●	●	●	●	●	●	●	●		
KÜSCHALL K-SERIES 2.0		●	●	●	-	-	-	-	-	-	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	●	●	●	●	-	●	●	●	●	●	●	●		
KÜSCHALL THE KSL 2.0		●	●	●	-	-	-	-	-	-	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	●	●	-	-	-	-	-	-	-	-	-	-	-	
KÜSCHALL CHAMPION version 1.0		-	●	●	-	●	-	-	-	-	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	
KÜSCHALL COMPACT, version 1.0	Design after 2009	●	●	●	●	●	-	-	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	
KÜSCHALL COMPACT ATTRACT		●	●	●	●	●	-	-	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	●	●	●	●	●	●	●	●	
KÜSCHALL K-SERIES, version 1.0 (G3)		●	●	●	-	-	-	-	-	-	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	
KÜSCHALL K-SERIES ATTRACT		●	●	●	-	-	-	-	-	-	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
KÜSCHALL THE KSL, version 1.0		-	●	●	-	-	-	-	-	-	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
KÜSCHALL ULTRA-LIGHT	Design after 2009	●	●	●	●	●	●	-	●	●	●	●	●	●	●	●	●	●	●	●	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	

*only compatible with fixed legrests - NC with swing away legrests

Combination Agreement



● Agreement compliant with MDR2017/745

Invacare manual wheelchair		Including versions	Bodypoint					Motion Concepts													Vicair								Supportec				
			FT240 Ankle huggers	NB2xxx-B5-ASM-50 Nonpadded Pelvic Positioning Belt	HB200 & HB400 series pelvic positioning belts using B1 end fittings	SH295 Harness Pivot Fit Dynamic	EB200 series EvoFlex Hip Belts	SP100 series Calf Supports	Backrest Matrix ELITE	Backrest Matrix ELITE DEEP	Backrest Matrix ELITE E2 (E2S)	Backrest Matrix ELITE E2 DEEP (E2D)	Backrest Matrix ELITE TR	Backrest Matrix PB POSTURE Back	Backrest Matrix PB DEEP POSTURE Back	Backrest Matrix MX2	Backrest Matrix Mini	Headrest Matrix ELAN and adapter	Headrest Matrix Mini ELAN	Headrest SOFT PAD and adapter Matrix / Neckrest Matrix	Cushion Matrix PS Paediatric	Cushion Matrix PS HD	Cushion Matrix-Vi	Active O2 seat cushion Adjuster O2 seat	Adjuster O2 seat cushion	Vicair Vector O2	Vicair Multi-functional O2 Back	Vicair Flex Back	Vicair Vector Junior O2	Vector Liberty Back	Vector Anatomic Back	Evoshape Soft Seat/Back Upholstered Evoshape Soft Seat/Back Parted Upholstered Evoshape Soft Flex Seat/Back Upholstered	
ACTION 1R			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
ACTION 2NG			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
ACTION 2 TRANSIT LITE			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
ACTION 3JUNIOR		Evolutive	●	●	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	
		Light	●	-	-	-	-	-	●	●	●	●	-	-	-	-	-	-	-	-	●	-	-	-	●	●	●	-	-	-	●	●	-
ACTION 3NG		Comfort	●	-	-	-	-	-	●	●	●	●	-	-	-	-	-	-	-	●	-	-	●	●	●	●	●	-	-	-	●	●	-
		XT	●	-	-	-	-	-	●	●	●	●	-	-	-	-	-	-	-	●	-	-	●	●	●	●	●	-	-	-	●	●	-
ACTION 3NG: ROCKING CHAIR			●	-	-	-	-	-	●	●	●	●	-	-	-	-	-	-	-	●	-	-	●	-	-	●	●	-	-	-	●	●	-
ACTION 4NG		Comfort	●	-	-	-	-	-	●	●	●	●	-	-	-	-	-	-	-	●	-	-	●	-	-	●	●	-	-	-	●	●	-
ACTION 4NG: HEAVY DUTY VERSION		XLT	●	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	●	-	-	-	●	●	-
ACTION 5/ MYON HC			-	-	-	-	-	-	●	●	●	●	-	-	-	-	-	-	-	-	-	-	●	-	-	-	-	-	-	-	-	-	-
ACTION 5/ MYON HC: RIGID VERSION			-	-	-	-	-	-	●	●	●	●	-	-	-	-	-	-	-	-	-	-	●	-	-	-	-	-	-	-	-	-	-
ACTION AMPLA			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-
ESPRIT ACTION 4NG			-	-	-	-	-	-	●	●	●	●	-	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-	-	-	-
ESPRIT ACTION			-	-	-	-	-	-	●	●	●	●	-	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-	-	-	-
SPIREA 4 NG			-	-	-	-	-	-	●	●	●	●	●	●	●	-	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-
		Standard	●	-	●	●	-	-	●	●	●	●	●	●	●	-	-	-	-	●	-	-	-	-	-	●	●	●	-	-	-	-	●
REA AZALEA		Assist	●	-	●	●	-	-	●	●	●	●	●	●	●	-	-	-	-	●	-	-	-	-	-	●	●	●	-	-	-	-	●
		Minor	-	-	-	-	-	-	●	●	●	●	●	●	●	-	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	●
		Tall	●	-	●	●	-	-	●	●	●	●	●	●	●	-	-	-	-	●	-	-	-	-	-	-	●	●	-	-	-	-	●
REA AZALEA BASE			-	-	-	-	-	-	●	●	●	●	●	●	●	-	-	-	-	●	-	-	-	-	-	-	-	-	-	-	-	-	-
REA AZALEA MAX			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
REA CLEMATIS			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-
REA CLEMATIS PRO		E-tilt	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●
REA DAHLIA		30, 45	●	-	●	●	-	-	●	●	●	●	●	-	-	-	-	-	-	●	-	-	-	-	-	●	●	●	-	-	-	-	●
KÜSCHALL CHAMPION 2.0			-	-	-	-	●	●	●	●	●	●	-	-	-	●	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-
KÜSCHALL COMPACT 2.0			-	-	-	-	●	●	●	●	●	●	-	-	-	●	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-
KÜSCHALL K-SERIES 2.0			-	-	-	-	●	●	●	●	●	●	-	-	-	●	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-
KÜSCHALL THE KSL 2.0			-	-	-	-	●	●	-	-	-	-	-	-	-	●	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-
KÜSCHALL CHAMPION version 1.0			-	-	-	-	-	-	●	●	●	●	-	-	-	●	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-
KÜSCHALL COMPACT, version 1.0		Design after 2009	-	-	-	-	-	-	●	●	●	●	-	-	-	●	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-
KÜSCHALL COMPACT ATTRACT			-	-	-	-	-	-	●	●	●	●	-	-	-	●	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-
KÜSCHALL K-SERIES, version 1.0 (G3)			-	-	-	-	-	-	●	●	●	●	-	-	-	●	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-
KÜSCHALL K-SERIES ATTRACT			-	-	-	-	-	-	●	●	●	●	-	-	-	●	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-
KÜSCHALL THE KSL, version 1.0			-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	●	-	-	-	-	-	-	-
KÜSCHALL ULTRA-LIGHT		Design after 2009	-	-	-	-	-	-	●	●	●	●	-	-	-	●	-	-	-	-	-	-	●	-	-	●	-	-	-	-	-	-	-