



Certificate in relation to air transportation by passengers with walking disabilities

We hereby certify that the our product

e-motion M25

including the lithium-ion batteries fulfil the requirements of

IATA chapter 2.3.2.4

and therefore may be carried as checked baggage without restriction.

The lithium-ion batteries of the e-motion M25 meet the requirements of the tests in the UN Manual of Tests and Criteria" Part III subsection 38.3.

The rated energy of the batteries in each wheel (**156,95 Wh**) is below the IATA limit of 160 Wh.

The e-motion M25 can be set into a dedicated "flight mode". In this mode the system is completely deactivated. While in this mode the e-motion M25 can only be switched on by a coded procedure. Unintended switching on by luggage is not possible. Furthermore the Bluetooth device of the system is completely switched off.

We strongly recommend to contact the airline sufficiently in advance to request approval for transport.

This certificate does not apply to damaged or defective batteries.

Due to the frequent amendments of the respective regulations, which partly occur on an annual basis this certificate is valid until **31st December 2019**.

Albstadt, 06th of February 2019

Alber GmbH

A handwritten signature in black ink, appearing to read "H. Reiner".

Hansjörg Reiner
Managing Director