



Product Service

CERTIFICATE

No. Z1 042097 0069 Rev. 00

Holder of Certificate: Invacare Rea AB
Org. No. 556098-9211
Växjövägen 303
343 71 Diö
SWEDEN

Certification Mark:



Product: Medical Beds

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

Test report no.: 713164526

Valid until: 2024-12-03

Date, 2020-02-27

(Torsten Zimmer)

TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD
ZERTIFIKAT ◆ CERTIFICATE ◆ 認證證書 ◆ CERTIFICADO ◆ CERTIFICAT



Product Service

CERTIFICATE

No. Z1 042097 0069 Rev. 00

Model(s): NordBed: Essential, Ultra,Optimo 90, Optimo wide 105, Optimo wide 120

Type	NordBed Essential	NordBed Ultra	NordBed Optimo 90	NordBed Optimo wide 105	NordBed Optimo wide 120
Application environment	3, 4	3, 4	3, 4	3, 4	3, 4
Safe working load	205 kg	235 kg	235 kg	235 kg	285 kg
Maximum patient weight	170 kg	200 kg	200 kg	200 kg	250 kg
Range of height adjustment	28 - 82 cm	25 - 85 cm	32 - 81 cm	32 - 81 cm	32 - 81 cm
Mattress support platform	90 x 200 cm	90 x 200 cm	90 x 200 cm	105 x 200 cm	120 x 200 cm
Voltage supply	230 VAC, 50 Hz				
Max. current input	3.5 A	4 A			
Operating time	max. 10%, 2 min/18 min				
Type of applied part	Type B				
Protection class	II				
Degree of protection	IP66W				

Tested according to:

- IEC 60601-1:2005
- IEC 60601-1:2005/AMD1:2012
- IEC 60601-1-6:2010
- IEC 60601-1-6:2010/AMD1:2013
- IEC 62366-1:2015
- IEC 60601-1-11:2015
- IEC 60601-2-52:2009
- IEC 60601-2-52:2009/AMD1:2015
- IEC 60601-1-2:2014

Production Facility(ies): 042097

TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD TÜV SÜD
 ZERTIFIKAT ◆ CERTIFICATE ◆ 認證證書 ◆ CERTIFICADO ◆ CERTIFICAT