







RULLSTOLENS HÖJDINSTÄLLNINGAR







hög 24x1^{3/8}"

						
41,5	3°	13	4	N/A	7	125
44	3°	11	2	N/A	5	125
44	3°	11	3	N/A	N/A	140
44	3°	11	3	N/A	6	150
46,5	3°	9	N/A	3	3	125
46,5	3°	9	1	4	N/A	140
46,5	3°	9	1	4	4	150
46,5	3°	9	2	N/A	N/A	180
49	3°	7	N/A	1	1	125
49	3°	7	N/A	2	N/A	140
49	3°	7	N/A	2	2	150
49	3°	7	N/A	3	N/A	180
49*	3°	7	1	4	4	200
51,5	3°	5	N/A	1	N/A	140
51,5	3°	5	N/A	1	1	150
51,5	3°	5	N/A	2	N/A	180
51,5*	3°	5	N/A	3	3	200
54*	3°	3	N/A	1	1	200

N/A = Ej tillämplig

* = Inte möjligt med 90° benstöd







hög 24x1"

						
41,5	3°	12	4	N/A	7	125
44	3°	10	2	N/A	5	125
44	3°	10	3	N/A	N/A	140
44	3°	10	3	N/A	6	150
46,5	3°	8	N/A	3	3	125
46,5	3°	8	1	4	N/A	140
46,5	3°	8	1	4	4	150
46,5	3°	8	2	N/A	N/A	180
49	3°	6	N/A	1	1	125
49	3°	6	N/A	2	N/A	140
49	3°	6	N/A	2	2	150
49	3°	6	N/A	3	N/A	180
49*	3°	6	1	4	4	200
51,5	3°	4	N/A	1	N/A	140
51,5	3°	4	N/A	1	1	150
51,5	3°	4	N/A	2	N/A	180
51,5*	3°	4	N/A	3	3	200
54*	3°	2	N/A	1	1	200

N/A = Ej tillämplig

* = Inte möjligt med 90° benstöd







hög 22^{3/8"}

						
41,5	3°	11	4	N/A	7	125
44	3°	9	2	N/A	5	125
44	3°	9	3	N/A	N/A	140
44	3°	9	3	N/A	6	150
46,5	3°	7	N/A	3	3	125
46,5	3°	7	1	4	N/A	140
46,5	3°	7	1	4	4	150
46,5	3°	7	2	N/A	N/A	180
49	3°	5	N/A	1	1	125
49	3°	5	N/A	2	N/A	140
49	3°	5	N/A	2	2	150
49	3°	5	N/A	3	N/A	180
49*	3°	5	1	4	4	200
51,5	3°	3	N/A	1	N/A	140
51,5	3°	3	N/A	1	1	150
51,5	3°	3	N/A	2	N/A	180
51,5*	3°	3	N/A	3	3	200
54*	3°	1	N/A	1	1	200

N/A = Ej tillämplig

* = Inte möjligt med 90° benstöd







hög 22x1"

						
41,5	3°	10	4	N/A	7	125
44	3°	8	2	N/A	5	125
44	3°	8	3	N/A	N/A	140
44	3°	8	3	N/A	6	150
46,5	3°	6	N/A	3	3	125
46,5	3°	6	1	4	N/A	140
46,5	3°	6	1	4	4	150
46,5	3°	6	2	N/A	N/A	180
49	3°	4	N/A	1	1	125
49	3°	4	N/A	2	N/A	140
49	3°	4	N/A	2	2	150
49	3°	4	N/A	3	N/A	180
49*	3°	4	1	4	4	200
51,5	3°	2	N/A	1	N/A	140
51,5	3°	2	N/A	1	1	150
51,5	3°	2	N/A	2	N/A	180
51,5*	3°	2	N/A	3	3	200
54*	3°	N/A	N/A	N/A	N/A	200

N/A = Ej tillämplig

* = Inte möjligt med 90° benstöd







hög 20x1^{3/8}"

						
41,5	3°	9	4	N/A	7	125
44	3°	7	2	N/A	5	125
44	3°	7	3	N/A	N/A	140
44	3°	7	3	N/A	6	150
46,5	3°	5	N/A	3	3	125
46,5	3°	5	1	4	N/A	140
46,5	3°	5	1	4	4	150
46,5	3°	5	2	N/A	N/A	180
49	3°	3	N/A	1	1	125
49	3°	3	N/A	2	N/A	140
49	3°	3	N/A	2	2	150
49	3°	3	N/A	3	N/A	180
49	3°	3	1	4	4	200

N/A = Ej tillämplig

* = Inte möjligt med 90° benstöd







Iåg 20x1^{3/8}"

						
36,5	3°	9	4	N/A	7	125
39	3°	7	2	N/A	5	125
39	3°	7	3	N/A	N/A	140
39	3°	7	3	N/A	6	150
41,5	3°	5	N/A	3	3	125
41,5	3°	5	1	4	N/A	140
41,5	3°	5	1	4	4	150
44	3°	3	N/A	1	1	125
44	3°	3	N/A	2	N/A	140
44	3°	3	N/A	2	2	150
44	3°	3	N/A	3	N/A	180

N/A = Ej tillämplig

* = Inte möjligt med 90° benstöd







Iåg 22x1^{3/8}"

						
36,5	3°	11	4	N/A	7	125
39	3°	9	2	N/A	5	125
39	3°	9	3	N/A	N/A	140
39	3°	9	3	N/A	6	150
41,5	3°	7	N/A	3	3	125
41,5	3°	7	1	4	N/A	140
41,5	3°	7	1	4	4	150
41,5	3°	7	2	N/A	N/A	180
44	3°	5	N/A	1	1	125
44	3°	5	N/A	2	N/A	140
44	3°	5	N/A	2	2	150
44	3°	5	N/A	3	N/A	180
44*	3°	5	1	4	4	200
46,5	3°	3	N/A	1	N/A	140
46,5	3°	3	N/A	1	1	150
46,5	3°	3	N/A	2	N/A	180
46,5*	3°	3	N/A	3	3	200
49*	3°	1	N/A	1	1	200

N/A = Ej tillämplig

* = Inte möjligt med 90° benstöd







låg 22x1"

						
36,5	3°	10	4	N/A	7	125
39	3°	8	2	N/A	5	125
39	3°	8	3	N/A	N/A	140
39	3°	8	3	N/A	6	150
41,5	3°	6	N/A	3	3	125
41,5	3°	6	1	4	N/A	140
41,5	3°	6	1	4	4	150
41,5	3°	6	2	N/A	N/A	180
44	3°	4	N/A	1	1	125
44	3°	4	N/A	2	N/A	140
44	3°	4	N/A	2	2	150
44	3°	4	N/A	3	N/A	180
44*	3°	4	1	4	4	200
46,5	3°	2	N/A	1	N/A	140
46,5	3°	2	N/A	1	1	150
46,5	3°	2	N/A	2	N/A	180
46,5*	3°	2	N/A	3	3	200
49*	3°	N/A	N/A	N/A	N/A	200

N/A = Ej tillämplig

* = Inte möjligt med 90° benstöd







låg 24x1^{3/8}"

						
39	3°	11	2	N/A	5	125
39	3°	11	3	N/A	N/A	140
39	3°	11	3	N/A	6	150
41,5	3°	9	N/A	3	3	125
41,5	3°	9	1	4	N/A	140
41,5	3°	9	1	4	4	150
41,5	3°	9	2	N/A	N/A	180
44	3°	7	N/A	1	1	125
44	3°	7	N/A	2	N/A	140
44	3°	7	N/A	2	2	150
44	3°	7	N/A	3	N/A	180
44*	3°	7	1	4	4	200
46,5	3°	5	N/A	1	N/A	140
46,5	3°	5	N/A	1	1	150
46,5	3°	5	N/A	2	N/A	180
46,5*	3°	5	N/A	3	3	200
49*	3°	3	N/A	1	1	200
49*	3°	N/A	N/A	N/A	N/A	200

N/A = Ej tillämplig

* = Inte möjligt med 90° benstöd

låg 24x1"

						
39	3°	10	2	N/A	5	125
39	3°	10	3	N/A	N/A	140
39	3°	10	3	N/A	6	150
41,5	3°	8	N/A	3	3	125
41,5	3°	8	1	4	N/A	140
41,5	3°	8	1	4	4	150
41,5	3°	8	2	N/A	N/A	180
44	3°	6	N/A	1	1	125
44	3°	6	N/A	2	N/A	140
44	3°	6	N/A	2	2	150
44	3°	6	N/A	3	N/A	180
44*	3°	6	1	4	4	200
46,5	3°	4	N/A	1	N/A	140
46,5	3°	4	N/A	1	1	150
46,5	3°	4	N/A	2	N/A	180
46,5*	3°	4	N/A	3	3	200
49*	3°	2	N/A	1	1	200

N/A = Ej tillämplig

* = Inte möjligt med 90° benstöd