

## COMBINATION AGREEMENT

Invacare REA AB  
Växjö 303  
343 71 DIÖ  
SWEDEN  
Orgnr: Orgnr:556618-9006

Decon Wheel AB  
Södra Ekeryd 115  
314 32 Hyltebruk  
Sweden

as manufacturers of the

Wheelchairs:

REA Spirea  
REA Focus  
REA SpinX  
REA Dahlia  
REA Clematis  
REA Azalea  
XLT  
REA Bellis

Accessories  
E-move  
E-drive  
E-walk  
Easy-Wheel

hereby declare that the above-mentioned products can be combined under the following conditions:

- The necessary fittings must be fitted onto the basic frame of the wheelchair.*  
Should it prove necessary to place the fittings on other parts of the wheelchair, the specific locations must be determined individually and agreed upon by both parties.
- Constructional changes to the wheelchair are allowed under following conditions.*  
Should it prove necessary to add or remove parts, or to make any other changes to the bed frame, such modifications must be approved of by the wheelchair manufacturer.
- The assembly instructions in the product manuals must be observed.*  
If any special fittings are required for specific wheelchair, the design and manufacturing thereof must be agreed upon by both parties. The mutually approved fittings must be accompanied by specific assembly instructions.
- Technical reservations concerning the combination must be specified.*  
Should either of the parties be of the opinion, that a combination requires specific equipment to be fitted onto a bed (e.g. wheels with certain dimensions) this must be agreed upon by both parties and stated in a specific assembly instruction.
- The unit must not limit any functions or adjustments on the wheelchair.*  
Certain changes to the initial adjustment possibilities of the bed are however acceptable, provided they have no influence on safety issues.
- Any further reservations in relation to the combination must be agreed upon by both parties.*

The units mentioned are CE-marked and meet the requirements of Council Directive 93/42/EEC concerning medical devices of Class I. Provided that the above conditions are met, the manufacturers hereby declare that the combination of the products also meets these requirements.

If second-hand products are to be used for the combination, they must first be inspected in accordance with their instructions for recycling. This also applies, if a unit is to be fitted onto one of the above-mentioned wheelchairs, which due to its age is not CE-marked.

If regulatory requirements should change or if constructional changes should be made to the products, each of the parties shall evaluate whether or not this affects the validity of the agreement.


Each of the parties to the agreement shall notify the other party of any such changes.

Should an incident occur with the combination, which would require notification of the national authorities, a copy of the agreement shall be enclosed the notification.

The agreement shall be valid for an indefinite period of time as from its execution. However, it can be terminated by either party in writing by giving 30 days notice.

9/12-2015

Date

  
Göran Löfgren


14/12-2015

Date

  
Benny Andersson

9/12-2015

Date

  
Gunnar Karlsson