



	DIMENSION RYGGSKAL UTAN SKUM			DIMENSION RYGGSKAL INCLUSIVE SKUM (KONTAKTYTA)				
	A BREDD	B HÖJD	C LATERAL KONTUR HÖJD	D INRE MAXIMAL BREDD	E YTTRE MAXIMAL BREDD	F HÖJD	VIKT UTAN MONTERINGSFÄSTE	MAX BRUKARVIKT
ARTIKELNUMMER	cm	cm	cm	cm	cm	cm	kg.	kg.
MINI1007	20	15	6	18	25	18	0.41	70
MINI1010	20	23	6	18	25	25	0.45	70
MINI1012	20	28	6	18	25	30	0.59	70
MINI1107	23	15	6	20	28	18	0.50	70
MINI1110	23	23	6	20	28	25	0.54	70
MINI1112	23	28	6	20	28	30	0.64	70
MINI1207	25	15	6	23	30	18	0.50	70
MINI1210	25	23	6	23	30	25	0.54	70
MINI1212	25	28	6	23	30	30	0.73	70
MINI1214	25	33	6	23	30	36	0.82	70
MINI1307	28	15	6	25	33	18	0.54	70
MINI1310	28	23	6	25	33	25	0.64	70
MINI1312	28	28	6	25	33	30	0.77	70
MINI1314	28	33	6	25	33	36	0.91	70
MINI1407	30	15	6	28	36	18	0.59	70
MINI1410	30	23	6	28	36	25	0.73	70
MINI1412	30	28	6	28	36	30	0.82	70
MINI1414	30	33	6	28	36	36	0.95	70